

BIOplastics

DeBox™ Family

Ticks all the Boxes

Unique work/storage system with endless possibilities!

- Sustainable
- Stackable
- Spotless



Get rid of cardboard boxes!

DeBox Family of BIOplastics consists of state-of-the-art stackable boxes and modular inserts to make any PCR-rack you like. This great work/storage solution ticks all the boxes and offers you:

176 boxes

Different dimensions, heights and hinge sides

Endless storage possibilities

Small boxes fit into medium size boxes and medium size boxes fit into large size boxes

Eight different colors

Create at least 792 different types of boxes

A wide variety of DeBox inserts

- Made from recycled and reprocessed Medical Grade PP
- Connectable to each other
- Indestructible

Products are designed and manufactured in The Netherlands / DeBox™ is a trademark of BIOplastics

© 2019 BIOplastics B.V. All rights reserved / 280519

DeBox Boxes


Mix & Match

Footprint	Dimensions	Available heights					
Small	85 x 59 mm	21.4 mm	35.2 mm	49 mm			
Medium	135 x 93 mm	26.7 mm	40.05 mm	53.4 mm			
Large	220.3 x 157 mm	32 mm	46 mm	60 mm	73 mm	87 mm	114 mm


Heights marked in the same color fit into each other

All boxes are available in:

- Natural, red, blue, green, yellow and white or a combination of two colors
- Hinges on the long or short side



Hinges long side



Hinges short side


Box in a box

Example 1

2 small boxes (85 x 59 x 21.4 mm) fit into
1 medium box (135 x 93 x 26.7 mm)
2 medium boxes (135 x 93 x 26.7 mm) fit into
1 large box (220.3 x 157 x 32 mm)




2 small boxes fit into a medium box



Example 2

2 small boxes (85 x 59 x 35.2 mm) fit into
1 medium box (135 x 93 x 40.05 mm)
2 medium boxes (135 x 93 x 40.05 mm) fit into
1 large box (220.3 x 157 x 46 mm)



2 medium boxes fit into a large box

Example 3

2 small boxes (85 x 59 x 49 mm) fit into
1 medium box (135 x 93 x 53.4 mm)
2 medium boxes (135 x 93 x 53.4 mm) fit into
1 large box (220.3 x 157 x 60 mm)



Products are designed and manufactured in The Netherlands / DeBox™ is a trademark of BIOplastics







© 2019 BIOplastics B.V. All rights reserved / 280519

DeBox Boxes







Mix & Match

Possible combinations

Small boxes (85 x 59 mm)

	Height	Hinges long side	Hinges short side
Low bottom Low lid	21.4 mm		
Medium bottom Low lid	35.2 mm		
Medium bottom Medium lid	49 mm		

Medium boxes (135 x 93 mm)

	Height	Hinges long side	Hinges short side
Low bottom Low lid	26.7 mm		
Medium bottom Low lid	40.05 mm		
Medium bottom Medium lid	53.4 mm		

Designed and manufactured in The Netherlands / DeBox™ is a trademark of BIOplastics













© 2019 BIOplastics B.V. All rights reserved / 280519

DeBox Boxes

Mix & Match

Possible combinations

Large boxes (220.3 x 157 mm)

	Height	Hinge long side	Hinge short side
Low bottom Low lid	32 mm		
Medium bottom Low lid	46 mm		
Medium bottom Medium lid	60 mm		
Large bottom Low lid	73 mm		
Large bottom Medium lid	87 mm		
Large bottom Large lid	114 mm		

Designed and manufactured in The Netherlands / DeBox™ is a trademark of BIOplastics

© 2019 BIOplastics B.V. All rights reserved / 280519


DeBox Inserts

Mix & Match

All DeBox inserts are connectable to each other and indestructible. You're no longer stuck to standard grids, but are able to create any platform you like, over and over again.


DeBox-In (q)PCR Insert 6 x 8 low

Dimensions in mm 55 x 81 x 10 (L x W x H)
 Capability (q)PCR vials, 8-tube strips and 48-well plate
 Colors natural, red, blue, green, yellow, white
 Package size box of 5
 Article **B11209** (white)
 Equals BLEGO 123 **BLEGO123-L9** (see next page)




DeBox-In (q)PCR Insert 1.5 - 4 x 8 | 32 x 0.2 ml + 3 x 1.5/2 ml + 2 x 0.5 ml

Dimensions in mm 55 x 81 x 20 (L x W x H)
 Capability (q)PCR vials, 8-tube strips and 32-well plate
 +3 x 1.5/2 ml & 2 x 0.5 ml tubes
 Colors natural, red, blue, green, yellow, white
 Package size box of 5
 Article **B11219** (white)




DeBox-In (q)PCR Insert 14 x 8 low

Dimensions in mm 128 x 86 x 10 (L x W x H)
 Capability (q)PCR vials, 8- and 12-tube strips & 96-well plate
 Colors natural, red, blue, green, yellow, white
 Package size box of 5
 Article **B11259** (white)



DeBox-In (q)PCR Insert 12 x 8 (0.2 ml) + 3 x 1.5/2 ml + 2 x 0.5 ml

Dimensions in mm 55 x 81 x 20 (L x W x H)
 Capability (q)PCR vials, 8- and 12-tube strips and 96-well plate
 +3 x 1.5/2 ml & 2 x 0.5 ml tubes
 Colors natural, red, blue, green, yellow, white
 Package size box of 5
 Article **B11269** (white)



Products are designed and manufactured in The Netherlands / DeBox™ is a trademark of BIOplastics


© 2019 BIOplastics B.V. All rights reserved / 280519

DeBox Inserts

Mix & Match


DeBox-In (q)PCR Insert Fantasy low

2 circles + 5 squares + 1 rectangle
 Dimensions in mm 55 x 81 x 10 (L x W x H)
 Colors white, natural
 Package size box of 5
 Article **B11229** (white)



DeBox-In (q)PCR Insert Fantasy high


2 circles + 5 squares + 1 rectangle
 Dimensions in mm 55 x 81 x 20 (L x W x H)
 Colors white, natural
 Package size box of 5
 Article **B11239** (white)



DeBox-In 123 has three interlocking elements. The blocks can be extended multiple times on all sides with any of the DeBox-inserts


BHLEGO 123 base DeBox-In (q)PCR 3 x 4

Dimensions in mm 72 x 18 x 24 (L x W x H) per element
 Capability 12 x 0.5, 1.5 or 2.0 ml microcentrifuge or screw cap tubes
 Color white, red, blue, yellow, green, natural
 Package size Pack of 7 (21 pieces)
 Article **BHLEGO123-L9** (white)



BLEGO 123 base DeBox-In (q)PCR 3 x 2 x 8

Dimensions in mm 72 x 18 x 10 (L x W x H) per element
 Capability 0.1 or 0.2 ml / 48-well (q)PCR plate or 6 x 8-tubes strips
 Color white, red, blue, yellow, green, natural
 Package size pack of 7 (21 pieces)
 Article **BLEGO123-L9** (white)



Equals DeBox-In (q)PCR Insert 6 x 8 low B11209 (see previous page)

BIOplastics

DeBox Family

Personalize...



Personalize your boxes and inserts!

Laser mark them with your logo, QR- or barcode, kit content,

You name it
We realize it!



Think out of the box
& create any rack you like



Products are designed and manufactured in The Netherlands/ DeBox™ is a trademark of BIOplastics

© 2019 BIOplastics B.V. All rights reserved / 280519

www.bioplastics.com

+31(0)455338750

order@bioplastics.com