

Cycler to Product Chart **PLATES** and Strips with Attached Caps

Cycler	Strips with Attached Caps				Plates											Plates 384 Format	
	CF78201	B79201	B72811	B72911	B70501	AB70651	B60101	B70651	B70671	AB17500	AB19800	B71501	B50501	B50651	B17489	B70515L	B71519L
	RP	HP	LP	RP	RP	LP	LP	RP	LP	RP	LP	RP	RP	RP	LP		
ABI																	
2400	V	V	X	V	Vc	X	X	X	X	X	X	Vc	Vc	Vc	X		
2700	V	V	X	V	V	X	X	V	X	V+	X	V	V	V	X	X	X
2720	V	V	X	V	V	X	X	V	X	V+	X	V	V	V	X	X	X
5700	V	V	X	V	V	X	X	V	X	V	X	V	V	V	X	X	X
7000	V	V	X	V	V	X	X	V	X	V	X	V	V	V	X	X	X
7300	V	X	X	V	V	X	X	X	X	V+	X	V	V+	V+	X	X	X
7500	V	X	X	V	V	X	X	X	X	V+	X	V	V+	V+	X	X	X
7500 Fast	X	X	Va	X	X	X	X	X	X	X	V+	X	X	X	X	X	X
7700	V	X	X	V	V	X	X	X	X	V+	X	V	V	V	X	X	X
7900	X	X	X	Va	X	X	X	X	X	V	X	V	V	V	X	V+	X
7900 Fast	X	X	Va	X	X	X	X	X	X	X	V+	X	X	X	X	X	X
9600	V	V	V	V	V	X	X	V	X	V+	X	V	V+	V+	V	X	X
9700	V	X	X	V	V	X	X	V	X	V+	X	V	V+	V+	X	V+	X
9800 Fast	X	X	V	X	X	V	V	X	X	X	V+	X	X	X	V	X	X
StepOne	X	X	V	X	X	Vc	Vc	X	X	X	X	X	X	X	Vc	X	X
StepOne Plus	X	X	V	X	X	V	V	X	X	X	V+	X	X	X	V	X	X
Veriti	V	X	X	V	V	X	X	V	X	V+	X	V	V+	V+	X	V+	X
Veriti Fast	X	X	V	X	X	V	V	X	X	X	V+	X	X	X	X	X	X
Biometra																	
Tpersonal	V	V	V	V	Vc	Vc	X	X	X	X	X	V	Vc	Vc	X	X	X
T1	V	V	V	V	V+	V	V	V	V	X	X	V	V	V	V	V	V
Tgradient	V	V	V	V	V+	V	V	V	V	X	X	V	V	V	X	X	X
T3000	V	V	V	V	Vc	X	X	X	X	X	X	V	Vc	Vc	X	X	X
Tprofessional	V	V	V	V	V+	V	V	V	V	X	X	V	V	V	X	V	V
Trobot	V	V	V	V	V+	V	V	V	V	X	X	V	V	V	V	V	V
BioRad																	
GeneCycler	V	V	X	V	Vc	X	X	V	X	X	X	Vc	Vc	Vc	X	X	X
iCycler	V	V	X	V	V	X	X	V	X	X	X	V	V+	V+	X	V	V
iCycler IQ	V	V	X	V	V	X	X	V	X	X	X	V	V+	V+	X	V	V
IQ5	V	V	X	V	V	X	X	V	X	X	X	V	V+	V+	X	X	X
MyCycler	V	V	X	V	V	X	X	V	X	X	X	V	V+	V+	X	X	X
MyIQ	V	V	X	V	V	X	X	V	X	X	X	V	V+	V+	X	X	X
Ptc100	V	V	V	V	V	V	V	V	V	X	X	V	V	V	V	V	V
Ptc150	X	X	V	V	X	X	X	X	X	X	X	X	X	X	X	X	X
Ptc200	V	V	V	V	V	V	V	V	V	X	X	V	V	V	V	V	V
Ptc225	V	V	V	V	V	V	V	V	V	X	X	V	V	V	V	V	V
Opticon	X	X	V	X	X	V	V	V	V	X	X	X	X	X	V	X	X
Opticon2	X	X	V	X	X	V	V	X	V	X	X	X	X	X	V	X	X
MiniOpticon	V	V	V	V	Vc	X	Vc	X	X	X	X	Vc	Vc	X	X	X	X
Chromo4	V	V	V	V	V	V	V	V	V	X	X	V	V	V	V	X	X
C-1000	X	X	V	X	X	V	V	X	V	X	X	X	X	V	V	X	X
F-1000	X	X	V	X	X	V	V	X	V	X	X	X	X	V	V	X	X
CFX96	V	X	V	V	V	V	V	V	V	X	X	V	V	V	V	X	X
CFX384	X	X	X	X	X	X	X	X	X	X	X	X	X	X	X	V	V
Eppendorf																	
MasterCycler	V	V	V	V	V	V	V	V	V	X	X	V	V	V	V	X	X
MC Gradient	V	V	V	V	V	V	V	V	V	X	X	V	V	V	V	X	X
MC ep Gradient	V	V	V	V	V	V	V	V	V	X	X	V	V	V	X	X	X
MC ep Gradient S	V	V	V	V	V	V	V	V	V	X	X	V	V	V	X	V	V
MC ep Realplex	V	V	V	V	V	V	V	V	V	X	X	V	V	V	X	V	V
MC Personal 16	X	X	X	X	X	X	X	X	X	X	X	X	X	X	X	X	X
Stratagene																	
Robo-Cycler	V	V	X	V	V	V	V	V	V	X	X	V	V	V	X	X	X
MX3000P	V	X	X	V	V	X	X	V	X	X	X	V	V	V	X	X	X
MX3005P	V	X	X	V	V	X	X	V	X	X	X	V	V	V	X	X	X
MX4000	X	X	X	V	V	X	X	V	X	X	X	V	V	V	X	X	X
Roche																	
LightCycler 480	X	X	Va	X	X	Va	Va	X	X	X	X	X	X	X	V+	X	V+
Lightlyper	X	X	Va	X	X	Va	Va	X	X	X	X	X	X	X	V+	X	V+
Cycler	RP	HP	LP	RP	RP	LP	LP	RP	LP	RP	LP	RP	RP	RP	LP		

Closures for qPCR (8 strips) 1, B57801 EU Wide Optical 8-Cap Strip
 2, B79701 EU Optical 8-Cap strip
 3, 157300 Optiseal
 Closures for PCR All strip caps and seals