




Compose your own (q)PCR plate by clicking the tube-strip into the grid

EU Thin Wall (q)PCR Robust 8-Tube Strip For Diagnostic & Clinical procedures

excellent performance and handling

- Perfect fit to nearly ALL PCR & qPCR cyclers
(see reference next page)
- Fits also ABI 7300 and 7500 Real Time Cyclers
- Designed to perform qPCR and PCR reactions
- Extra robust, no twisting tubes or cycler placing or removal problems
- Low evaporation due to clever design
- Easy labeling and orientation
- Extreme Uniform tube leading to reproducible results
- Can be closed using any of the BIOplastics cap strips
- Can be placed, transported & stored in Tube support grid, station support insert or Work Rack S-96



EU 0.2ml Thin-wall 8-tube strip **Extra Robust**

Can be used in all 0.2 ml block regular cyclers and most Real-Time thermal cyclers (see table next page). Closure can be accomplished with EU optical, flat cap thin-wall 8-cap strips (B79701, B75701, B57801) or with the EU optical, semi-domed 8-cap strips (C79701). For Real time applications B57801 is recommended

<u>Orderno</u>	<u>Description</u>	<u>Package Size</u>	<u>Price</u>
B69901	EU 0.2ml Thin-wall 8-tube strip Extra Robust	Natural, Bag of 120 Strips	
B57801	EU Wide Optical 8-Cap Strip, Wide Indented Flat Cap	Natural, Bag of 120 Strips	
B75701	EU Optical 8-Cap strip	Natural, Bag of 120 Strips	
B69351	EU Pre-Post Tube Support Grid Wide, Natural, Re-useable	Natural, 8 Grids	
B69303	EU Pre-Post Tube Support Grid, Blue, Re-useable	Blue, 8 Grids	



q)PCR Cycler reference chart for Regular Profile robust 8-tube strip


All Cycler brands with height adjustable lids: YES



	310	3100	3130	3700	3730	2400	9600	2700	2720	7000	7300	7500	7700	7900	9700	Veriti™ 96
Applied Biosystems (non fast)	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓
	GeneCycler	iCycler	MyCycler	PTC 100	PTC 200	PTC 220	PTC 221	PTC 225	Chromo-4	iCycler IQ	IQ 5	MyiQ	CFX 96	C1000	S1000	
Biorad	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓
	MasterCycl.	Master Grad.	Master ep	ep Realplex												
Eppendorf	✓	✓	✓	✓												
	Mx 3000P	Mx 3005P														
Stratagene	✓	✓														



B72711 (natural) & B72719 (white)																	
Low Profile robust strips which fits																	
	7500 fast*	7900 fast*	9800	StepOne	StepOne Plus	Verity Fast	Opticon I	Opticon II	Mini Opticon	Roche LC 480*							
																	* Specific adaptor required



Products are designed and manufactured in the Netherlands. All rights reserved. 012009V1
 Trademarks; owner: marks- > Applied Biosystems: - ABI Prism - Verity - StepOne - StepOne Plus ; Bio-Rad: - Opticon - iCycler - MyCycler - Chromo-4 - C1000 - S1000 ; Eppendorf: - MasterCycler - RealPlex ep ; Stratagene: - Mx3000P - Mx3005P

Head Quarters:
 (design and Manufacturing)
 BIOplastics BV
 Röttscherweg 61
 Entrance D
 6374 XW Landgraaf
 The Netherlands

 Phone: +31 (0)45 533 87 50
 Fax: +31 (0)45 533 87 96
 Website: www.bioplastics.com
 E-mail: info@bioplastics.com

Subsidiary:

 BPCTi Inc.
 2933 South Miami Blvd.
 Suite 121
 Durham, NC 27703, USA

 Phone: 919 806 8811
 Fax: 919 806 2014
 E-mail: info@bpcti.com
 Website: www.bpcti.com

Your local distributor: