

Work rack systems

Get OrgaNICEd

Work | Store | Freeze

- Autoclavable
- Stackable
- Durable



Characteristics

BIOplastics work rack systems are the perfect solution to work with, store and freeze microcentrifuges tubes & cryo vials. Temperature range from $-180\text{ }^{\circ}\text{C}$ to $110\text{ }^{\circ}\text{C}$. They are made from durable polypropylene so extremely strong. The boxes and racks are stackable so you can organise your lab as you wish.

Hygienic

The boxes and work racks are autoclavable & can not be saturated with water like cardboard boxes.

Content examination

It is possible to view the content of the box without removing the lid.

Colours

The boxes and some of the work racks are available in up to eight different colors. In this product info sheet the article numbers refer to blue boxes with natural lid and white work racks.

Packaging

All articles are packed per eight boxes/work racks/work rack systems, unless mentioned otherwise.

Designed and manufactured in The Netherlands. Samples on request.

© 2019 Bioplastics B.V. All rights reserved / 100719

Work rack systems

Work rack S System



Box S regular model


Dimensions box in cm: 13.1 (L) x 9.2 (W) x 6.2 (H)

Also available in 10.5 cm high model: [B10133](#)

[B10113](#)


Fitting work racks

Dimensions in cm: 12.3 (L) x 8.4 (8.4(W) x 4.5 (H)




[B10079](#)

Work rack S-40
For 40 x 1.5 ml & 2 ml microcentrifuge tubes



[B10099](#)


Work rack S-60
For 60 x 0.5 ml microcentrifuge tubes



[B69409](#)


Work rack S-96
For 96 x 0.2 ml tubes, strips or plates

Combinations



[B10053](#)

Box & work rack S-40



[B10033](#)

Box & work rack S-60



[B10013](#)


Box & work rack S-96

Designed and manufactured in The Netherlands. Samples on request.

© 2019 Bioplastics B.V. All rights reserved / 100719

Work rack systems

Work rack F100 System




Box F100

Dimensions box in cm: 14.3 (L) x 14.3 (W) x 5.9 (H)

B10313

Fitting work racks


Dimensions in cm: 13.2 (L) x 13.2 (W) x 4.4 (H)



B10259

Work rack F81


For 81 x 1.5 ml & 2 ml microcentrifuge tubes



B10279

Work rack F combo

96 x 0.2 ml tubes/strips/plates & 16 x 1.5/2 ml tubes & 8 x 0.5 ml tubes




B10239

Work rack F100


100 x 0.5 ml tubes

Combinations




B10173

Box & work rack F-81



B10193

Box & work rack F combo



B10153


Box & work rack F100

Designed and manufactured in The Netherlands. Samples on request.

© 2019 Bioplastics B.V. All rights reserved / 100719

Work rack systems

Multo work rack system



Multo work rack box
B10403

Dimensions box in cm:
10.36 (L) x 13.5 (W) x 2.77 (H)



Fitting work rack
B10429

Multo plate rack holds low- & regular profile tubes, strips & plates.
Laser marked A-H & 1-12




Combination
B10443

Box & work rack Multo

Medical kits

The Multo box closes very securely. Due to its limited height of less than 3 cm (1.2 inch) it's a perfect solution for packaging medical kits or transporting valuable samples.



Cryo storage box 81

B10021

The cryo box has 81 positions & a working range from ambient temperatures from -200 °C to 110 °C
Packed per 5 boxes

Modular work/storage solutions

For modular work/storage solutions, you can check out the newest **DeBox Family** of BIOplastics



Designed and manufactured in The Netherlands. Samples on request.

© 2019 Bioplastics B.V. All rights reserved / 100719